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

Smart Heart Rate and Blood Pressure Wristband

User Guide


Quick instructions for use

Part of introduction

Charging method

For the first time, please take the wristband out and charge the equipment with the corresponding charger. When charging, the wristband will open automatically. (charge mode as follows: Joint the USB-charging cable with Band at match point backsideand USB port for charging.)

Notice: according to the style of the bracelet, please select the above figure A or B for disassembly

Quick instructions for use

Power on/off

1. In the state of the wristband shutdown, the long touch function key over 3 seconds can make the device vibrate.
2. When the wristband is in the boot state, it will switch to the shutdown interface. If the touch function key is more than 3 seconds, it can enter the shutdown interface. After the short touch is selected Off, it will touch the shutdown again.

operation

1. When the device starts up, to touch function keys for a short time period can lighten screen or switch over display;
2. The screen shall be turned off without operation and users can connect the badminton interface show calories, long press three seconds to exit the sport. mobile phone APP to revise off-screen time;
3. Switch to the heart rate and blood pressure interfaces to be tested automatically and leave the heart rate and blood pressure interfaces to stop tests.

Install the wristband APP


Scan the following QR codes or enter into various large application markets to download and install APP

Device requirement: iOS8.0 and above; Android 4.4 and above to support bluetooth 4.0 .

Device attachment

When being used for the first time, it is necessary to connect APP for alignment. After successful connection, wristband will self-synchronize time, or otherwise, step calculation and sleep data won't be accurate.

Bluetooth connection: after successfully matching, APP will automatically save the wristband blueooth address. When APP is opened or operated in background, it will automatically search and link with the wristband.

Data Synchronism: manually pull down the synchronization data on the APP home page; the bracelet can store off-line data for 7 days, the more data, the longer the synchronization time is, and the longest time is about 2 minutes.

Main functional interface

Main interface

display real-time temperature, bluetooth connection indicator icon, power indicator icon, charging indicator picture, time, date, week, percentage of movement steps.

Stepping interface

Display step count, wristband default step counting mode, the wristband is worn on the wrist. When the user walks and arms swings, the wristband is recorded according to the number of steps.

Distance interface

The interface is based on the number of walking steps, the height of the APP individual, and the journey of the day.

Calorie interface

The interface is based on the number of walking steps, the height and weight of the APP individual, and calculates the total number of calories consumed on the day.

Heart rate interface

At the heart rate monitoring interface, it automatically enters the heart rate detection mode and supports dynamic real-time heart rate detection. In the heart rate test, the duration of the wristband will be reduced.

Blood pressure interface

Switch to the blood pressure interface, start testing blood pressure, switch to the next interface to stop testing, start testing, but do not measure the result, show the results of last test or " ---/.. " ; also connect to mobile phone

App, and test at the blood pressure interface. Test methods and attention items of concentric rate detection.

Blood oxygen interface

Switch to the blood oxygen pressure interface, start testing blood oxygen pressure, switch to the next interface to stop testing, start testing, but do not measure the result, show the results of last test or " ---/.. " ; also connect to mobile phone App, and test at the blood pressure interface. Test methods and attention items of concentric rate detection.

Sports interface

Walking mode: long press to start walking movement:
Running mode: start running with long press:
In walking and running mode, the page is displayed by two pages The first page
1. Step 2. The distance 3. The calories 4. A stopwatch
The second page
1. The clock 2. Average speed 3. The heart rate 4. A stopwatch
Data display page long press 1S to end less movement mode and display motion summary information.

Message interface

open the message history, according to historical information entry, can save up to 16 the history of the latest news, long press 1 seconds can enter to check the news page, and then through the short press recycled switch to the next message, in any YIYITiao history news page long press 1 seconds to quit, the last page to delete all the trash mutations news page, long press 1 seconds to delete all message history.

Find phone interface (optional)

When connected, click on the "Find Phone" option and the phone will send a vibration reminder or ringtone reminder.

Version information

displays the version number and MAC address

Shutdown interface

long press more than three seconds to shake the phone off

APP functions and setting

- 1) APP reminder: WeChat, QQ, phone, SMS, facebook messenger, twitter reminder to open/close control and reserve other message alerts.
- 2) alarm reminder: 5 sets of smart alarm clocks can be set, which can be ordinary alarm clock (display the alarm time) or alarm clock (display the information entered by the user). Repeat alarm and single alarm are supported.
- 3) sedentary reminder: turn on the sedentary reminder, which can set the time and duration of sitting, and trigger the reminder beyond the set time.
- 4) do not disturb mode: turn on the do not disturb mode and you can set the do not disturb time period. During this time period, except the alarm will vibrate, all other reminders will be closed.
- 5) look for bracelet: when connected, look for bracelet. The bracelet displays the look up icon and vibrates.
- 6) hand lifting screen: support the function of hand lifting screen after opening.
- 7) measurement cycle: heart rate automatically measures the time cycle. This setting is effective when the measurement switch is turned on.
- 8) sleep time: you can set the sleep time, during which the alarm will vibrate except for the alarm clock. All other alarms will be closed. In addition, the sleep time and sleep quality of human

body will be recorded during this period. During sleep time, if the heart rate monitoring is not turned on, the sleep is monitored during the whole sleep time. If the heart rate monitoring is turned on, the user is first monitored to see if he is wearing it. If he is wearing it, he will start to monitor his sleep. If he is not wearing it, he will not record his sleep.

- 9) firmware upgrade: OTA wireless upgrade is supported. If the firmware version of the APP query server bracelet is higher than the current version, firmware upgrade can be selected.
- 10) firmware version: current version of the bracelet
- 11) bluetooth address: bracelet MAC address
- 12) disconnection reminder: it is suggested to open the broken bracelet after opening to prevent loss.
- 13) take photos of the bracelet: Enter the bracelet to take photos when connected with bluetooth. The bracelet displays the photo page.
- 14) exercise mode: turn on the screen to display the movement information directly, and turn off the screen to display the first main page when the screen is turned off. Turn on the bracelet to display the screen and automatically measure the heart rate every 15 minutes.
- 15) restore factory Settings: clear all RAM and FLASH information and data, and reset the bracelet.
- 16) restart the bracelet: clear RAM data and reset the bracelet.
- 17) store the history of the bracelet: when opened, you can store the received message notice, up to 16 messages can be stored, and the earliest and latest messages can be discarded after more than 16 messages. Turning this feature off does not save the message.
- 18) delete the bracelet history message: the whole history message of the bracelet can be deleted by APP operation (the bracelet can be deleted by pressing "delete all history message" on the last page of the history message).

Real-time synchronization of bracelet data: the number of steps, distance and calories of the bracelet can be synchronized in real time after being opened (real-time data is displayed on the APP terminal).

My feature set

- 1) set nickname; Sex; Age; Height; Weight. The above information can be used to calculate the walking distance and burn calories. When the bracelet is not set, the default is: male, 30 years old, 177cm, 80KG
- 2) my MEDALS and my weekly newspaper show statistics
- 3) target movement: the number of target steps can be set
- 4) setting: clear the mobile phone data (note: the bracelet can store valid data for 14 days, and the mobile phone can always save data if it is not clear)

running function

After opening, the track and distance of movement can be recorded

home page function

- 1) show daily steps, calories, exercise distance, target completion percentage, average heart rate, blood pressure, blood oxygen, and daily sleep time (deep sleep and light sleep time) .
- 2) click to see the bar chart of the number of steps
- 3) click to view the measured data of heart rate, blood pressure and blood oxygen. Single measurement (1 minute) and real-time measurement (10 minutes) are supported in the heart rate data page. At this time, the measured data will be displayed on the APP in real time.
- 4) click to view the distribution of sleep quality, time of deep sleep and number of awakenings.